



ChemDraw/Signals ChemDraw 23.0



Efficient Chemical Drawing



Chemical & Biological Intelligence



Eye Catching Visual Chemistry



Integrations with Key Research Applications & Customizations



Quick & Painless Chemistry Reporting



Modern Cloud Architecture

ChemDraw Prime



- Comprehensive Chemical Drawing Tools
- Hotkeys & Shortcuts
- Generic Structures
- Structure & Reaction Clean-up
- Smart Copy/Paste (SMILES to Structure, auto-positioning)



- Enhanced stereo perception

ChemDraw Professional



- Comprehensive Chemical Drawing Tools
- Hotkeys & Shortcuts
- Generic Structures
- Structure & Reaction Clean-up
- Smart Copy/Paste (SMILES to Structure, auto-positioning)



- ¹H & ¹³C NMR Prediction
- Name-to-Structure
- Structure-to-Name
- HELM Biopolymer Toolbar
- Centralized biopolymer library management in Signals (HELM Monomer Curation)
- Enhanced stereo perception



- 3D Clean-up
- Atom / Bond Color Highlighting
- Ring-Fill Coloring
- Structure Perspective Depth-fading
- Biopolymer Colored Document Settings



- SciFinder-n Integration
- Reaxys Integration
- Millipore Sigma Synthia Integration

Signals ChemDraw



- Comprehensive Chemical Drawing Tools
- Hotkeys & Shortcuts
- Generic Structures
- Structure & Reaction Clean-up
- Smart Copy/Paste (SMILES to Structure, auto-positioning)



- ¹H & ¹³C NMR Prediction
- Name-to-Structure
- Structure-to-Name
- HELM Biopolymer Toolbar
- Centralized biopolymer library management in Signals (HELM Monomer Curation)
- Enhanced stereo perception



- [Save and Copy as 3D Printable Format \(.3MF\)](#)
- [Transfer Color Highlight to 3D Printable Format \(.3MF\)](#)
- 3D Clean-up
- Atom / Bond Color Highlighting
- Ring-Fill Coloring
- Structure Perspective Depth-fading
- Biopolymer Colored Document Settings



- [Google Scholar/Patents Add-in](#)
- [PubChem Safety Add-in](#)
- [ChemACX Explorer Add-in](#)
- [ChemDraw JavaScript License for Internal Development](#)
- [Custom, Dynamic Add-in SDK](#)
- [Centralized HELM Monomer Libraries](#)
- SciFinder-n Integration
- Reaxys Integration
- Millipore Sigma Synthia Integration



- ChemDraw Collections
- Cloud Anywhere Access and Automatic Updates
- Create Collections of Molecules/Reactions
- Add/Edit Properties in Collections
- Export Collections to .sdf/.pptx Files
- Create Reports of Signals Notebook Experiments
- Extract ChemDraw/ISIS Draw Drawings from Documents
- Perform Structure Searches in Cloud-stored Documents



- Centralized user management (SaaS)
- ChemDraw Desktop automatic updates
- Create and manage & curate Chemical drawings in the cloud (ChemDraw+)

Revvity Signals Software, Inc
940 Winter Street | Waltham, MA 02451 USA
P: (800) 762-4000 or (+1) 203-925-4602
revvitysignals.com/company/contact

Copyright ©, Revvity, Inc. All rights reserved. Revvity™ is a registered trademark of Revvity, Inc. All other trademarks are the property of their respective owners.

revvity
signals