

InventoryGo: Mobile Inventory Management for Signals Notebook and Signals One

revvity
signals

By Revvity Signals

February 2026

Quick Facts

Audit-Ready Inventory Management

InventoryGo enables accurate, compliant, mobile-first inventory management within Signals One and Signals Notebook. It supports lab efficiency, traceability, and audit readiness while minimizing IT and security concerns.

Introduction

Fragmented, manual inventory systems create inefficiencies, stockouts, and compliance risks for scientific labs and procurement teams. Customers need a reliable, centralized inventory system that integrates with scientific workflows to reduce manual tracking, ensure accurate stock levels, support regulatory compliance, and allow hands-free use in glove-worn lab environments.

Revvity Signals InventoryGo is a lab inventory app that allows scientists to perform real-time inventory actions directly at the bench. Purpose-built for scientific environments, InventoryGo improves accuracy, compliance, and efficiency while integrating seamlessly with Signals Notebook and Signals One.

The Benefits of InventoryGo

InventoryGo is a progressive web app (PWA) that provides scientists with mobile-first access to inventory capabilities, including real-time reconciliation and inventory consumption during experiments, minimizing the need to operate computers with gloved hands. The first mobile extension of Signals Notebook and Signals One designed specifically for laboratory inventory management, InventoryGo installs on iOS, Android, Windows, and Mac devices directly—no app store required.

- Accurate inventory updates: Delivered in real-time at the bench.
- Reduced reliance on desktop computers: Ideal for gloved environments.
- Improved regulatory compliance: Strengthens traceability and audit readiness.
- Enhanced operational efficiency: No need to introduce new infrastructure or APIs.
- Purpose-built UI: Designed to be intuitive, low-maintenance, and impactful with minimal training.

Who is it for?

InventoryGo supports bench scientists, lab managers, and procurement teams who rely on Signals One and Signals Notebook and want fast, reliable access to inventory during active experiments.

Have Confidence in Your Lab Inventory

- Install anywhere, fast: PWA is available on iOS, Android, Windows, Mac.
- Reliable lab inventory: Realtime stock updates.
- Kiosk mode: Facilitates check-in and check-out using shared devices.
- Scientific workflow integration: Purpose-built for science, InventoryGo improves accuracy, compliance, and efficiency.
- Secure access: Leverage existing Signals One and Signals Notebook permissions.

InventoryGo: Built into Your Scientific Process

Deeply Integrated with the Signals One Solution


InventoryGo is fully integrated with Signals One and Signals Notebook, so users will enjoy consistent permissions, data governance, and compliance. It requires no additional infrastructure or new APIs.

Device-Agnostic Compatibility

InventoryGo is offered as a Progressive Web App (PWA), delivering a fast, app-like experience across mobile, tablet, and desktop devices. The PWA provides native-style speed, offline-ready behavior, and responsive performance across platforms. It can currently be accessed in two ways:

- **Through a web browser**, by logging into Signals One or Signals Notebook and launching the InventoryGo PWA, which loads directly within the browser and provides applike responsiveness.
- **Through a downloadable mobile application**, installed directly on supported devices for users who prefer a native app experience.

Together, these access options ensure broad compatibility and a seamless experience, regardless of device or operating system.



Revvity Signals' broader Inventory capabilities

Reliable and intuitive inventory management wherever science happens, the bench, storage rooms, or on the move. [Request a Demo Today!](#)

Purpose-Built for Scientific Environments

InventoryGo has a reliable, fast, and engaging UI designed for scientists. It offers the following inventory actions:

- **Container:** Check out or move containers.
- **Consume:** Update container amounts in real time.
- **Location:** Check containers in and out of specific storage locations.
- **Reconciliation:** Reconcile and rectify discrepancies instantly.

To learn more about Revvity Signals' broader Inventory capabilities, visit [Inventory | Revvity Signals Software](#)

revvity signals

77 4th Avenue
Waltham, MA 02451 USA
(800) 762-4000
(+1) 203-925-4602
revvitysignals.com



Copyright © Revvity, Inc. All rights reserved. Revvity™ is a registered trademark of Revvity, Inc. All other trademarks are the property of their respective owners.